



Endothelial Cell Medium-rat-basal-phenol red free (ECM-r-b-prf)

Catalog #1001-r-b-prf

Product Description

Endothelial Cell Medium-rat-basal-phenol red free (ECM-r-b-prf), when used with Endothelial Cell Growth Supplement-rat (ECGS-r, Cat #1062) and 25 ml of calf serum (CS, Cat #0025-cs) is a complete medium designed for optimal growth of normal rat endothelial cells *in vitro*. It is a sterile, liquid medium which contains essential and non-essential amino acids, vitamins, organic and inorganic compounds and trace minerals. The medium is HEPES and bicarbonate buffered and has a pH of 7.4 when equilibrated in an incubator with an atmosphere of 5% CO₂/95% air.

Components

ECM-r-b-prf consists of 500 ml of basal medium.

Product Use

ECM-r-b-prf is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

Storage

Store at 4°C. Protect from light.

Shipping

Room temperature.

Caution: If handled improperly, some components of the medium may present a health hazard. Take appropriate precautions when handling it, including the wearing of protective clothing and eyewear. Dispose of properly.