



Oligodendrocyte Precursor Cell Medium-basal (OPCM-b)

Catalog #1601-b

Product Description

Oligodendrocyte Precursor Cell Medium-basal (OPCM-b), when used with Oligodendrocyte Precursor Cell Growth Supplement (OPCGS, Cat #1652), is a complete medium designed for optimal growth of normal oligodendrocyte precursor cells *in vitro*. It is a sterile, liquid medium which contains essential and non-essential amino acids, vitamins, organic and inorganic compounds and trace minerals. The medium is HEPES and bicarbonate buffered and has a pH of 7.4 when equilibrated in an incubator with an atmosphere of 5% CO₂/95% air.

Components

OPCM-basal consists of 500 ml of basal medium.

Product Use

OPCM-b is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

Storage

Store the basal medium at 4°C and the OPCGS and P/S solution at -20°C. Protect from light.

Shipping

Room temperature.

may present a health hazard. Take appropriate precautions when handling it, including the wearing of protective clothing and eyewear. Dispose of properly.