



**Oligodendrocyte Precursor Cell Differentiation Medium-basal
(OPCDM-b)**
Catalog #1631-b

Product Description

Oligodendrocyte Precursor Cell Differentiation Medium-basal (OPCDM-b), when used with Oligodendrocyte Precursor Cell Differentiation Supplement (OPCDS, Cat #1672) and 5 ml of fetal bovine serum (FBS) is a complete medium designed for optimal growth of normal oligodendrocyte precursor cells *in vitro*. It is a sterile, liquid medium which contains essential and non-essential amino acids, vitamins, organic and inorganic compounds and trace minerals. The medium is HEPES and bicarbonate buffered and has a pH of 7.4 when equilibrated in an incubator with an atmosphere of 5% CO₂/95% air.

Components

OPCDM-b consists of 500 ml of basal medium.

Product Use

OPCDM-b is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

Storage

Store at 4°C. Protect from light.

Shipping

Room temperature.

Caution: If handled improperly, some components of the medium may present a health hazard. Take appropriate precautions when handling it, including the wearing of protective clothing and eyewear. Dispose of properly.