



## **Keratinocyte Medium-animal component free-basal (KM-acf-b)**

Catalog #2121-b

### **Product Description**

Keratinocyte Medium-animal component free-basal (KM-acf-b), when used with Keratinocyte Growth Supplement-animal component free (KGS-acf, Cat #2172), is a complete medium designed for optimal growth of normal keratinocytes *in vitro*. It is a sterile, liquid medium which contains essential and non-essential amino acids, vitamins, organic and inorganic compounds and trace minerals. The medium is HEPES and bicarbonate buffered and has a pH of 7.4 when equilibrated in an incubator with an atmosphere of 5% CO<sub>2</sub>/95% air.

### **Components**

KM-acf-basal consists of 500 ml of basal medium.

### **Product Use**

KM-acf-b is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

### **Storage**

Store the basal medium at 4°C and the KGS-acf and P/S solution at -20°C. Protect from light.

### **Shipping**

Room temperature.

*may present a health hazard. Take appropriate precautions when handling it, including the wearing of protective clothing and eyewear. Dispose of properly.*