



## **Fibroblast Medium-serum free-basal-phenol red free (FM-sf-b-prf)**

Catalog #2311-b-prf

### **Product Description**

Fibroblast Medium-serum free-basal-phenol red free (FM-sf-b-prf), when used with Fibroblast Growth Supplement-serum free (FGS-sf, Cat #2362), is a complete medium designed for optimal growth of normal fibroblasts *in vitro*. It is a sterile, liquid medium which contains essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors and trace minerals. The medium is HEPES and bicarbonate buffered and has a pH of 7.4 when equilibrated in an incubator with an atmosphere of 5% CO<sub>2</sub>/95% air.

### **Components**

FM-sf-b-prf consists of 500 ml of basal medium.

### **Product Use**

FM-sf-b-prf is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

### **Storage**

Store at 4°C. Protect from light.

### **Shipping**

Room temperature.

*Caution: If handled improperly, some components of the medium may present a health hazard. Take appropriate precautions when handling it, including the wearing of protective clothing and eyewear. Dispose of properly.*