



**Human Colonic Microvascular Endothelial Cell Pellet
(HCoMEC Cell Pellet)**
Catalog #CP2930

Quantity: 5 million cells

Product Description

Human Colonic Microvascular Endothelial Cell Pellet (HCoMEC Cell Pellet) is prepared from early passage Human Colonic Microvascular Endothelial Cells. Each pellet contains 5 million cells and can be used for a variety of applications including PCR, western blotting, genomic DNA library construction, and gene expression profiling. Cell pellets from ScienCell Research Laboratories are convenient and cost effective for researchers as it eliminates the need to acquire expensive tissues.

Applications

- Genomic DNA
- cDNA
- Total RNA
- microRNA
- Protein Lysates

Product Use

HCoMEC Cell Pellet is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

Storage

Store cell pellet immediately at -80°C upon arrival. Avoid freeze/thaw cycles.

Shipping

Dry ice.