



**Human Pulmonary Small Airway Epithelial Cell Pellet
(HPSAEpiC Cell Pellet)**
Catalog #CP3230

Quantity: 5 million cells

Product Description

Human Pulmonary Small Airway Epithelial Cell Pellet (HPSAEpiC Cell Pellet) is prepared from early passage Human Pulmonary Small Airway Epithelial Cell. Each pellet contains 5 million cells and can be used for a variety of applications including PCR, western blotting, genomic DNA library construction, and gene expression profiling. Cell pellets from ScienCell Research Laboratories are convenient and cost effective for researchers as it eliminates the need to acquire expensive tissues.

Applications

- Genomic DNA
- cDNA
- Total RNA
- microRNA
- Protein Lysates

Product Use

HPSAEpiC Cell Pellet is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

Storage

Store cell pellet immediately at -80°C upon arrival. Avoid freeze/thaw cycles.

Shipping

Dry ice.