



**Human Seminal Vesicle Microvascular Endothelial Cell Pellet  
(HSVMEC Cell Pellet)**  
Catalog #CP4450

Quantity: 5 million cells

**Product Description**

Human Seminal Vesicle Microvascular Endothelial Cell Pellet (HSVMEC Cell Pellet) is prepared from early passage Human Seminal Vesicle Microvascular Endothelial Cells. Each pellet contains 5 million cells and can be used for a variety of applications including PCR, western blotting, genomic DNA library construction, and gene expression profiling. Cell pellets from ScienCell Research Laboratories are convenient and cost effective for researchers as it eliminates the need to acquire expensive tissues.

**Applications**

- Genomic DNA
- cDNA
- Total RNA
- microRNA
- Protein Lysates

**Product Use**

HSVMEC Cell Pellet is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

**Storage**

Store cell pellet immediately at -80°C upon arrival. Avoid freeze/thaw cycles.

**Shipping**

Dry ice.