



**Human Cardiac Microvascular Endothelial Cell Pellet
(HCMEC Cell Pellet)**
Catalog #CP6000

Quantity: 5 million cells

Product Description

Human Cardiac Microvascular Endothelial Cell Pellet (HCMEC Cell Pellet) is prepared from early passage Human Cardiac Microvascular Endothelial Cells. Each pellet contains 5 million cells and can be used for a variety of applications including PCR, western blotting, genomic DNA library construction, and gene expression profiling. Cell pellets from ScienCell Research Laboratories are convenient and cost effective for researchers as it eliminates the need to acquire expensive tissues.

Applications

- Genomic DNA
- cDNA
- Total RNA
- microRNA
- Protein Lysates

Product Use

HCMEC Cell Pellet is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

Storage

Store cell pellet immediately at -80°C upon arrival. Avoid freeze/thaw cycles.

Shipping

Dry ice.