



**Human Cardiac Fibroblast-adult atrial Cell Pellet
(HCF-aa Cell Pellet)**
Catalog #CP6320

Quantity: 5 million cells

Product Description

Human Cardiac Fibroblast-adult atrial Cell Pellet (HCF-aa Cell Pellet) is prepared from early passage Human Cardiac Fibroblasts-adult atrial. Each pellet contains 5 million cells and can be used for a variety of applications including PCR, western blotting, genomic DNA library construction, and gene expression profiling. Cell pellets from ScienCell Research Laboratories are convenient and cost effective for researchers as it eliminates the need to acquire expensive tissues.

Applications

- Genomic DNA
- cDNA
- Total RNA
- microRNA
- Protein Lysates

Product Use

HCF-aa Cell Pellet is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

Storage

Store cell pellet immediately at -80°C upon arrival. Avoid freeze/thaw cycles.

Shipping

Dry ice.