



**Human Pericardial Fibroblast Cell Pellet  
(HPcF Cell Pellet)**  
Catalog #CP6430

Quantity: 5 million cells

**Product Description**

Human Pericardial Fibroblast Cell Pellet (HPcF Cell Pellet) is prepared from early passage Human Pericardial Fibroblasts. Each pellet contains 5 million cells and can be used for a variety of applications including PCR, western blotting, genomic DNA library construction, and gene expression profiling. Cell pellets from ScienCell Research Laboratories are convenient and cost effective for researchers as it eliminates the need to acquire expensive tissues.

**Applications**

- Genomic DNA
- cDNA
- Total RNA
- microRNA
- Protein Lysates

**Product Use**

HPcF Cell Pellet is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

**Storage**

Store cell pellet immediately at -80°C upon arrival. Avoid freeze/thaw cycles.

**Shipping**

Dry ice.