



Human Non-Pigmented Ciliary Epithelial Cell Pellet
(HNPCEpiC Cell Pellet)
Catalog #CP6580

Quantity: 5 million cells

Product Description

Human Non-Pigmented Ciliary Epithelial Cell Pellet (HNPCEpiC Cell Pellet) is prepared from early passage Human Non-Pigmented Ciliary Epithelial Cells. Each pellet contains 5 million cells and can be used for a variety of applications including PCR, western blotting, genomic DNA library construction, and gene expression profiling. Cell pellets from ScienCell Research Laboratories are convenient and cost effective for researchers as it eliminates the need to acquire expensive tissues.

Applications

- Genomic DNA
- cDNA
- Total RNA
- microRNA
- Protein Lysates

Product Use

HNPCEpiC Cell Pellet is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

Storage

Store cell pellet immediately at -80°C upon arrival. Avoid freeze/thaw cycles.

Shipping

Dry ice.